

## » Content

Implantswiss Dental System	2
Swiss Precision	4
Product Traceability	5
Glass Package	6
Platform Color Coding	7
<b>» Implantswiss System Specifications</b>	8
S.R.A. Surface	12
Hybrid Design	13
Conical Connection with Hexagonal Index	14
Platform Switching	14
Self Cutting Edge	15
Hygienic Neck Design	15
<b>» Implantswiss Implants</b>	16
Bone Level Implant	18
Tissue Level Implant	19
<b>» Surgical Protocols</b>	20
Bone Level-Tissue Level	21
Implantswiss Surgical Drills	22
<b>» Short Implant</b>	26
Short Implant	29
Short Implant Surgical Drills	30
<b>» Zygoma Implant</b>	36
Zygoma Implant Surgical Drills	37
<b>» Superb Overdenture System</b>	38
<b>» Implantswiss Prosthetic Components</b>	39
<b>» Bone Level System Prosthetic Components</b>	40
Bone Level Healing Abutments	42
Couple System	42
Solid System	44
Octa System	45
O-ring System	46
Locator System	47
Superb Overdenture System	48
Multi-Unit System	48
Bar System	49
<b>» Tissue Level System Prosthetic Components</b>	50
Tissue Level Healing Abutments	52
Couple System	52
Solid System	54
Octa System	55
O-ring System	56
Locator System	57
Superb Overdenture System	58
<b>» Zygoma Dental Implant System Prosthetic Components</b>	59
<b>» Drivers</b>	60
<b>» Screws</b>	61
<b>» Code Index</b>	62

implantswiss

Switzerland

is at the heart of

# science, mechanical design and craftsmanship



The well respected and acknowledged world authority...

Switzerland

has been one of the most important centers for craftsmanship that requires

# high precision and expertise

throughout history.

The products with impeccable mechanical designs, can be used for a lifetime and add value to your life.

**Implantswiss Dental System** owes its **excellent** structure to the inspiration, experience and value from history.

**Implantswiss** is a strong and reliable brand in the dental industry with its innovative solutions, unique surface design and high quality.

Implantswiss is devoted to Research and Development through cooperations with universities and crucial investments in clinical research. Implantswiss follows the latest trends and initiates extensive and contemporary projects that may lead the future of implant dentistry.

# Implantswiss Dental System


is the most advanced stage  
in implant dentistry



# SWISS PRECISION




Implantswiss products are manufactured with the state-of-the-arts technologies and high precision. Implantswiss never compromises on its quality and reliability; the success and high quality of Implantswiss lies in this perfectionist approach in each step of the manufacturing process.



# PRODUCT TRACEABILITY

Each Implantswiss product can be easily traced from production to delivery.

Each product package has a barcode label with a lot number created by the ERP software system. All steps from production to stock, delivery and after sales can be monitored with the barcode system; therefore, the follow-up process can be conducted flawlessly.



# GLASS PACKAGE

## GLASS PACKAGE

The nature and human health have always been the primary concern for Implantswiss. Implantswiss, gives great importance to protection of nature and human health. Therefore, all Implantswiss Implants are packed in 100% natural and hygienic glass vials.

Since the glass structure doesn't contain any foreign material, glass packages are especially preferred for preserving the high value products.

The high barrier property, ease of sterilisation and unreactivity in chemical substances are among the many advantages of glass packages.

## SMART PACKAGE


The implant is simply removed from the glass bottle and placed in the osteotomy using the fixture driver.

Each product package has a label containing all necessary information such as the type, the length and the color code of the product.

The color coding and associated diameters of the implants are classified as given.



# PLATFORM COLOR CODING




TISSUE LEVEL

BONE LEVEL



» **IMPLANTSWISS  
SYSTEM  
SPECIFICATIONS**






# » SRA\* Surface


Following calcium phosphate blasting, the surface of the implants is etched using Double Acid Etching.

Implantswiss S.R.A. surface has a homogenous surface roughness that perfectly mimics the bone morphology.

\* Sandblasted Rough Acid-etched Surface




## 3D-Area Parameters



Detail of the thread Bank and representation of the surface shape after removal of a 4th order polynomial

ISO 25178				
Height Parameters				
Sa	1,80	µm		
Sq	2,44	µm		
Feature Parameters				
S10z	15,1	µm	pruning = 5%	Ten point height
S5p	8,59	µm	pruning = 5%	Five point peak height
S5v	6,49	µm	pruning = 5%	Five point of height




Detail of the thread valley and representation of the surface shape after removal of a 4th order polynomial

ISO 25178				
Height Parameters				
Sa	1,25	µm		
Sq	1,60	µm		
Feature Parameters				
S10z	9,16	µm	pruning = 5%	Ten point height
S5p	3,39	µm	pruning = 5%	Five point peak height
S5v	5,76	µm	pruning = 5%	Five point of height

» Tapered form + Straight form = **Hybrid Design**

It is a unique implant design with cylindrical (straight) design on the neck area and conical (tapered) design through the apex of the implant.



The first 3 mm of Implantswiss implants running straight downward from the neck constitutes the straight form. The remaining portion constitutes the tapered form at a certain angle, depending on the length of the implant. The optimum combination of these two specifications creates the hybrid form. One advantage of the hybrid


design is to exert minimal stress on crestal cortical bone because of the straight portion at the neck, leading to an optimal load absorption rate. Due to its tapered form, 1/3 of the implant length is directly positioned in the osteotomy, which facilitate easy insertion of the implant in place.

## » Conical connection with hexagonal index

creates an excellent implant-abutment connection.

Implantswiss implants create an excellent implant-abutment interface, which prevents the formation of micro-gap and micro-leakage thanks to the universal conical connection and hexagonal index.

In Implantswiss Dental Implant System, all abutments are compatible with all fixtures with any diameter (except for Q 3.3) making a user friendly implant system for dentists.



## » Platform Switching

protects the alveolar bone around the neck area.

Micro-gap can be formed at the implant-abutment interface in Bone Level implants especially if they have a butt joint connection. Bacteria and by-products in these areas can cause serious problems such as bone resorption and periimplantitis.

Even though Implantswiss Implants have bacteria sealed conical connection, Switch Platform system further protects cervical bone by moving the implant-abutment interface away from alveolar bone. Implantswiss Implants with Switch Platform System are safer against bacteria and microleakage and have a longer life span.



## » Self Cutting Edge

Implantswiss implants have a thread design which does not transfer stress to the bone during insertion. The cutting edges and self cutting feature provide ideal primary stability in any bone type without transferring stress thanks to its semi-aggressive thread design.

Unique thread design that does not transfer stress to the bone.




## » Hygienic Neck Design

There are no micro threads on the neck area of Implantswiss implants therefore it is easier to maintain the hygiene when required.

Bacteria can easily colonize in the micro threads. In case of supra-alveolar positioning of the implant neck, the bone resorption risk is therefore minimized.

Easy to maintain the hygiene



» **IMPLANTSWISS  
IMPLANTS**




# BONE LEVEL IMPLANT



Bone Level implants are designed to provide excellent outcome in any bone quality and in any location.

Implantswiss bone level implants allow to design an ideal smile in the aesthetic zone.




8.0 mm	S-BMFSR3308	S-BFHR3708	S-BFHR4308	S-BFHR4808	S-BWFSR5508
10.0 mm	S-BMFSR3310	S-BFHR3710	S-BFHR4310	S-BFHR4810	S-BWFSR5510
12.0 mm	S-BMFSR3312	S-BFHR3712	S-BFHR4312	S-BFHR4812	S-BWFSR5512
14.0 mm	S-BMFSR3314	S-BFHR3714	S-BFHR4314	S-BFHR4814	-



# TISSUE LEVEL IMPLANT


An excellent soft and hard tissue protection can be achieved with the smooth neck and location of implant-abutment interface of the Tissue Level implant.

The implant is designed to provide excellent outcomes in any case where the aesthetic outcome is not the primary concern.



## » SURGICAL PROTOCOLS

BONE LEVEL - TISSUE LEVEL




# IMPLANTSWISS FIXTURE DRILLS




## SURGICAL DRILLS

Self stoppered surgical drills eliminate the risk of complications both during the preparation of the osteotomy and insertion of the implant. Implantswiss drills are recommended to be used at 200-300 rpm to prevent over heating and reduce mechanical trauma, which will shorten the osseointegration time of the implant. The drills are designed to harvest autogenous bone when used below 50rpm without irrigation to be used as a graft material in regenerative procedures.

Implantswiss Surgical Kit contains a specific drill for each implant available in the system.




Implantswiss Surgical Kit includes various drills with different diameters and lengths. These drills are designed to simplify the surgical procedures resulting in shortened surgery time and decreased stress for the patients.




**ø3.3  
Fixture**




D1 Bone



D2-3 Bone




D4 Bone




**ø3.7  
Fixture**




**D1 Bone**



**D2-3 Bone**




**D4 Bone**




**ø4.3  
Fixture**




**D1 Bone**



**D2-3 Bone**




**D4 Bone**




**ø4.8  
Fixture**




**D1 Bone**



**D2-3 Bone**



**D4 Bone**





**ø5.5  
Fixture**



**D1 Bone**

Point drill	2.2 drill	Pilot drill	3.3 Fixture drill	3.7 Fixture drill	4.3 Fixture drill	4.8 Fixture drill	5.5 Fixture drill	5.5 Neck drill	5.5 Fixture
-------------	-----------	-------------	-------------------	-------------------	-------------------	-------------------	-------------------	----------------	-------------

**D2-3 Bone**

Point drill	2.2 drill	Pilot drill	3.3 Fixture drill	3.7 Fixture drill	4.3 Fixture drill	4.8 Fixture drill	5.5 Fixture drill	5.5 Fixture
-------------	-----------	-------------	-------------------	-------------------	-------------------	-------------------	-------------------	-------------

**D4 Bone**

Point drill	2.2 drill	Pilot drill	3.3 Fixture drill	3.7 Fixture drill	4.3 Fixture drill	4.8 Fixture drill	5.5 Fixture
-------------	-----------	-------------	-------------------	-------------------	-------------------	-------------------	-------------

» **IMPLANTSWISS  
SHORT IMPLANTS**




# SHORT IMPLANT


In cases of insufficient bone volume, either above the mandibular canal or below the maxillary sinuses, short implants may be an alternative treatment to advanced surgical interventions aiming to regenerate hard tissues.

Due to aggressive thread design, Bone level short implants can be preferred mainly in maxilla, where the alveolar bone is relatively soft, whereas Tissue level short implant can be a good choice especially for the lower jaw.

## BONE LEVEL



## TISSUE LEVEL




# SHORT IMPLANT SURGICAL PROTOCOL

The surgical kit is recommended to be used with both Bone and Tissue Level short implants.



# IMPLANTSWISS SHORT IMPLANT DRILLS

Implantswiss Short Implant Surgical Kit contains a specific drill for each implant available in the system.



Implantswiss Short Implant Surgical Kit includes various drills with different diameters and lengths. These drills are designed to simplify the surgical procedures resulting in shortened surgery time and decreased stress for the patients.




**ø4.3  
Fixture**




**BONE LEVEL**

D4 / D3 Bone



**TISSUE LEVEL**

D2 / D1 Bone




BONE LEVEL SHORT IMPLANTS ARE recommended to be used only in D3/D4 and Tissue level implants only in D1/D2 bone.

**ø4.8  
Fixture**




**BONE LEVEL**

D4 / D3 Bone



**TISSUE LEVEL**

D2 / D1 Bone






**ø5.5  
Fixture**




**BONE LEVEL**

D4 / D3 Bone



**TISSUE LEVEL**

D2 / D1 Bone



» **IMPLANTSWISS  
ZYGOMA IMPLANTS**



# ZYGOMA IMPLANT

Zygomatic implants are designed for cases where traditional implants can not be used and for patients who want to benefit from immediate loading protocols even in the presence of minimum residual maxillary bone height.

Implantswiss Zygomatic Implant offers three different product options.

## ZYGOMA IMPLANT

ø4.3



35 mm	S-ZYG35
37.5 mm	S-ZYG37.5
40 mm	S-ZYG40
42.5 mm	S-ZYG42.5
45 mm	S-ZYG45
47.5 mm	S-ZYG47.5
50 mm	S-ZYG50
52.5 mm	S-ZYG52.5
55 mm	S-ZYG55
57.5 mm	S-ZYG57.5
60 mm	S-ZYG60

## ATROPHIC MAXILLARY ZYGOMA IMPLANT

ø4.3



35 mm	S-ZAM35
37.5 mm	S-ZAM37.5
40 mm	S-ZAM40
42.5 mm	S-ZAM42.5
45 mm	S-ZAM45
47.5 mm	S-ZAM47.5
50 mm	S-ZAM50
52.5 mm	S-ZAM52.5
55 mm	S-ZAM55
57.5 mm	S-ZAM57.5
60 mm	S-ZAM60


## ONCOLOGICAL ZYGOMA IMPLANT

ø4.3



25 mm	S-ONC25
27.5 mm	S-ONC27.5
30 mm	S-ONC30
32.5 mm	S-ONC32.5

# ZYGOMA IMPLANT SURGICAL PROTOCOL




## SURGICAL DRILLS

Point drill	Pilot drill 3.5	Pilot drill Long	Twist drill 2.5 Short	Twist drill 2.5 Long	Twist drill 2.9 Short	Twist drill 2.9 Long	Twist drill 3.5 Short	Twist drill 3.5 Long
-------------	--------------------	---------------------	--------------------------	-------------------------	--------------------------	-------------------------	--------------------------	-------------------------



# SUPERB OVERDENTURE SYSTEM



## Functionality

Structure for distributing vertical and lateral pressure. (specially designed with combination of abutment and housing)

## Aesthetics

Excellent aesthetics achieved without discoloration of acrylic from conventional metal housing.

## Reliability

Contact between the head of the abutment with the ceramic balls within the housing with minimal abrasion.
















## Flexibility

Allows up to 40 degrees angulation between implants without loss of retention force.















» **IMPLANTSWISS  
PROSTHETIC COMPONENTS**

# BONE LEVEL

## COUPLE SYSTEM
















FLOW	FIXTURE	HEALING	ABUTMENT	LAB ANALOG	ABUTMENT	SCREW
 <p>ø3.3    ø3.7    ø4.3</p>  <p>ø4.8    ø5.5</p>	 <p>Closing Screw</p>  <p>Healing Abutment</p>	<p>Pick-up</p>  <p>Hex    N-Hex</p> <p>Transfer</p>  <p>Hex    N-Hex</p>		<p>Couple</p>  <p>Hex    N-Hex</p> <p>Multi</p>  <p>Hex    N-Hex</p> <p>Angled</p>  <p>15°    25°</p>		
TOOL	 <p>Ratchet Driver</p>	 <p>Hex Screw Driver</p>	 <p>Hex Screw Driver</p>			 <p>Hex Screw Driver</p>

## SOLID SYSTEM











FLOW	FIXTURE	HEALING	ABUTMENT	IMPRESSION CAP	LAB ANALOG	PLASTIC COPING
 <p>ø3.3    ø3.7    ø4.3</p>  <p>ø4.8    ø5.5</p>	 <p>Closing Screw</p>  <p>Healing Abutment</p>	<p>Solid</p>  <p>Solid</p> <p>Cap</p>  <p>Cap</p>			<p>Single</p>  <p>Single</p> <p>Bridge</p>  <p>Bridge</p>	
TOOL	 <p>Ratchet Driver</p>	 <p>Hex Screw Driver</p>	 <p>Solid Driver</p>  <p>Hex Screw Driver</p>			 <p>Hex Screw Driver</p>



# OCTA SYSTEM





FLOW	FIXTURE	HEALING	SUB OCTA ABUTMENT	IMPRESSION COPING	LAB ANALOG	LAB CYLINDER	SCREW
 <p>ø3.3    ø3.7    ø4.3</p>  <p>ø4.8    ø5.5</p>	 <p>Closing Screw</p>  <p>Healing Abutment</p>	 <p>Octa Healing Cap</p>  <p>Octa Abutment</p>	 <p>Pick-up</p>  <p>Transfer</p>		 <p>Octa</p> <p>N-Octa</p>		
TOOL	 <p>Ratchet Driver</p>	 <p>Hex Screw Driver</p>	 <p>Octa Driver</p>				 <p>Hex Screw Driver</p>

# O-RING SYSTEM

FLOW	FIXTURE	HEALING	O-RING ABUTMENT	O-RING LAB ANALOG	RETAINER
 <p>ø3.3    ø3.7    ø4.3</p>  <p>ø4.8    ø5.5</p>	 <p>Closing Screw</p>  <p>Healing Abutment</p>			 <p>O-Ring</p> <p>O-Ring Retainer</p>	
TOOL	 <p>Ratchet Driver</p>	 <p>Hex Screw Driver</p>	 <p>O-Ring Driver</p>		

# BONE LEVEL HEALING ABUTMENTS




## BONE LEVEL HEALING ABUTMENTS

				
	ø3.7*	ø4.5	ø5.5	ø6.5
<b>1 mm</b>	-	S-BHA4510	S-BHA5510	S-BHA6510
<b>2 mm</b>	-	S-BHA4520	S-BHA5520	S-BHA6520
<b>2.5 mm</b>	S-BMHA3725	-	-	-
<b>3 mm</b>	-	S-BHA4530	S-BHA5530	S-BHA6530
<b>3.5 mm</b>	S-BMHA3735	-	-	-
<b>4 mm</b>	-	S-BHA4540	S-BHA5540	S-BHA6540
<b>5 mm</b>	S-BMHA3750	S-BHA4550	S-BHA5550	S-BHA6550
<b>7 mm</b>	S-BMHA3770	-	-	-

\*only for ø3.3 implant

**System** » Bone Level  
**Material** » Ti Grade 5 ELI  
**Driver** » Hex Screw Driver  
**Torque Value** » 15-20 Ncm





## BONE LEVEL SHORT IMPLANT HEALING ABUTMENTS

			
	ø5.0	ø6.0	ø8.0
<b>2 mm</b>	S-BMHA5020	S-BMHA6020	S-BMHA8020
<b>3 mm</b>	S-BMHA5030	S-BMHA6030	S-BMHA8030

**System** » Bone Level  
**Material** » Ti Grade 5 ELI  
**Driver** » Hex Screw Driver  
**Torque Value** » 15-20 Ncm

# BONE LEVEL COUPLE SYSTEM




## BONE LEVEL COUPLE ABUTMENTS

				
	ø3.7*	ø4.5	ø5.5	ø6.5
<b>1 mm</b>	S-BMCA37107	S-BCA45107	S-BCA55107	S-BCA65107
<b>2 mm</b>	S-BMCA37207	S-BCA45207	S-BCA55207	S-BCA65207
<b>3 mm</b>	S-BMCA37307	S-BCA45307	S-BCA55307	S-BCA65307
<b>4 mm</b>	S-BMCA37407	S-BCA45407	S-BCA55407	S-BCA65407
<b>5 mm</b>	S-BMCA37507	S-BCA45507	S-BCA55507	S-BCA65507

\*only for ø3.3 implant





**System** » Bone Level  
**Material** » Ti Grade 5 ELI  
**Torque Value** » 25 Ncm  
**Driver** » Hex Screw Driver

## BONE LEVEL SHORT IMPLANT COUPLE ABUTMENTS

			
	ø5.0	ø6.0	ø8.0
<b>2 mm</b>	S-BMCA5020	S-BMCA6020	S-BMCA8020
<b>3 mm</b>	S-BMCA5030	S-BMCA6030	S-BMCA8030

**System** » Bone Level  
**Material** » Ti Grade 5 ELI  
**Torque Value** » 25 Ncm  
**Driver** » Hex Screw Driver





## BONE LEVEL COUPLE ABUTMENTS NON-HEX

	 ø3.7*	 ø4.5	 ø5.5	 ø6.5
<b>1 mm</b>	S-BMCAN37107	S-BCAN45107	S-BCAN55107	S-BCAN65107
<b>2 mm</b>	S-BMCAN37207	S-BCAN45207	S-BCAN55207	S-BCAN65207
<b>3 mm</b>	S-BMCAN37307	S-BCAN45307	S-BCAN55307	S-BCAN65307
<b>4 mm</b>	S-BMCAN37407	S-BCAN45407	S-BCAN55407	S-BCAN65407
<b>5 mm</b>	S-BMCAN37507	S-BCAN45507	S-BCAN55507	S-BCAN65507

\*only for ø3.3 implant

**System** » Bone Level  
**Material** » Ti Grade 5 ELI  
**Torque Value** » 25 Ncm  
**Driver** » Hex Screw Driver





## BONE LEVEL ANGLED ABUTMENT

	 ø3.7*	 ø4.5 15° 25°	 ø5.5 15° 25°	 ø6.5 15° 25°
<b>1 mm (15°)</b>	S-BMAA37151	<b>2 mm 15°</b> S-BAA45152	S-BAA55152	S-BAA65152
<b>2 mm (15°)</b>	S-BMAA37152	<b>2 mm 25°</b> S-BAA45252	S-BAA55252	S-BAA65252
<b>4 mm (15°)</b>	S-BMAA37154	<b>4 mm 15°</b> S-BAA45154	S-BAA55154	S-BAA65154
		<b>4 mm 25°</b> S-BAA45254	S-BAA55254	S-BAA65254

\*only for ø3.3 implant

**System** » Bone Level // **Material** » Ti Grade 5 ELI // **Torque Value** » 25 Ncm // **Driver** » Hex Screw Driver

## BONE LEVEL ANGLED ABUTMENT NON-HEX

	 ø3.7*	 ø4.5 15° 25°	 ø5.5 15° 25°	 ø6.5 15° 25°
<b>1 mm (15°)</b>	S-BMAAN37151	<b>2 mm 15°</b> S-BAAN45152	S-BAAN55152	S-BAAN65152
<b>2 mm (15°)</b>	S-BMAAN37152	<b>2 mm 25°</b> S-BAAN45252	S-BAAN55252	S-BAAN65252
<b>4 mm (15°)</b>	S-BMAAN37154	<b>4 mm 15°</b> S-BAAN45154	S-BAAN55154	S-BAAN65154
		<b>4 mm 25°</b> S-BAAN45254	S-BAAN55254	S-BAAN65254

\*only for ø3.3 implant

**System** » Bone Level // **Material** » Ti Grade 5 ELI // **Torque Value** » 25 Ncm // **Driver** » Hex Screw Driver

## BONE LEVEL MULTI ABUTMENTS // HEX / NON-HEX

	 ø3.7	 ø4.5	 ø5.5
<b>HEX</b>	S-BMMA37	S-BMA45	S-BMA55
<b>NON-HEX</b>	S-BMMAN37	S-BMAN45	S-BMAN55

**System** » Bone Level  
**Material** » Ti Grade 5 ELI  
**Torque Value** » 25 Ncm  
**Driver** » Hex Screw Driver


## TEMPORARY ABUTMENTS

$\varnothing 3.7^*$      $\varnothing 4.5$      $\varnothing 5.5$   
 S-BMTA37    S-BTA45    S-BTA55

\*only for  $\varnothing 3.3$  implant

**System** » Bone Level  
**Material** » PEEK  
**Torque Value** » 20 Ncm  
**Driver** » Hex Screw Driver

## BONE LEVEL IMPRESSION COPING PICK UP



	$\varnothing 3.7^*$	$\varnothing 4.5$	$\varnothing 5.5$	$\varnothing 6.5$
<b>HEX</b>	S-BMICP37L S-BMICP37S	S-BICP45L S-BICP45S	S-BICP55L S-BICP55S	S-BICP65L S-BICP65S
<b>NON-HEX</b>	S-BMICPN37L S-BMICPN37S	S-BICPN45L S-BICPN45S	S-BICPN55L S-BICPN55S	S-BICPN65L S-BICPN65S

\*only for  $\varnothing 3.3$  implant

**System** » Bone Level  
**Material** » Ti Grade 5 ELI

## BONE LEVEL IMPRESSION COPING TRANSFER



	$\varnothing 3.7^*$	$\varnothing 4.5$	$\varnothing 5.5$	$\varnothing 6.5$
<b>HEX</b>	S-BMICT37L S-BMICT37S	S-BICT45L S-BICT45S	S-BICT55L S-BICT55S	S-BICT65L S-BICT65S
<b>NON-HEX</b>	S-BMICTN37L S-BMICTN37S	S-BICTN45L S-BICTN45S	S-BICTN55L S-BICTN55S	S-BICTN65L S-BICTN65S

\*only for  $\varnothing 3.3$  implant

**System** » Bone Level  
**Material** » Ti Grade 5 ELI

## BONE LEVEL FIXTURE LAB ANALOG




S-BMFLA    S-BFLA

only for  $\varnothing 3.3$  implant

**System** » Bone Level  
**Material** » Ti Grade 5 ELI

# BONE LEVEL SOLID SYSTEM

## BONE LEVEL SOLID ABUTMENTS



	$\varnothing 3.7^*$	$\varnothing 4.5$	$\varnothing 5.5$	$\varnothing 6.5$
<b>1 mm</b>	S-BMSA3710	S-BSA4510	S-BSA5510	S-BSA6510
<b>2 mm</b>	S-BMSA3720	S-BSA4520	S-BSA5520	S-BSA6520
<b>3 mm</b>	S-BMSA3730	S-BSA4530	S-BSA5530	S-BSA6530
<b>4 mm</b>	S-BMSA3740	S-BSA4540	S-BSA5540	S-BSA6540
<b>5 mm</b>	S-BMSA3750	S-BSA4550	S-BSA5550	S-BSA6550

\*only for  $\varnothing 3.3$  implant

**System** » Bone Level  
**Material** » Ti Grade 5 ELI  
**Torque Value** » 25 Ncm  
**Driver** » Hex Screw Driver

## BONE LEVEL SOLID ABUTMENT CAP



ø3.7\*

S-BMSACap37

\*only for ø3.3 implant



ø5.1

S-BSACap45



ø5.8

S-BSACap55



ø6.8

S-BSACap65

**System** » Bone Level Solid System  
**Material** » POM

## BONE LEVEL SOLID ABUTMENT IMPRESSION CAP



ø3.7\*

S-BMSAICap37

\*only for ø3.3 implant



ø4.5

S-BSAICap45



ø5.5

S-BSAICap55



ø6.5

S-BSAICap65

**System** » Bone Level Solid System  
**Material** » Plastic

## BONE LEVEL SOLID LAB ANALOG



ø3.7\*

S-BMSALA37

\*only for ø3.3 implant



ø4.5

S-BSALA45



ø5.5

S-BSALA55



ø6.5

S-BSALA65

**System** » Bone Level Solid System  
**Material** » Ti Grade 5 ELI

# BONE LEVEL OCTA SYSTEM

## BONE LEVEL OCTA ABUTMENTS



ø3.7\*

**1 mm**

S-BMOA10



ø4.8

**2 mm**

S-BMOA20

S-BOA20

**3 mm**

S-BMOA30

S-BOA30

**4 mm**

S-BMOA40

S-BOA40

**5 mm**

S-BMOA50

S-BOA50

\*only for ø3.3 implant

**System** » Bone Level  
**Material** » Ti Grade 5 ELI  
**Torque Value** » : Abutment Screw  
25 Ncm Oclusal Screw 20 Ncm  
**Driver** » Octa Driver

## BONE LEVEL OCTA ABUTMENT HEALING CAP



ø4.8

**4 mm**

S-OAHCap

**6 mm**

S-OAHCapL

**System** » Octa System  
**Material** » Ti Grade 5 ELI

## BONE LEVEL OCTA SYSTEM PLASTIC CYLINDER



**OCTA** S-OAPICy  
**NON-OCTA** S-OAPICyN

**System** » Octa System  
**Material** » POM

## BONE LEVEL OCTA SYSTEM TITANIUM CYLINDER

**HEX** S-OATICy  
**NON-HEX** S-OATICyN

**System** » Octa System  
**Material** » Ti Grade 5 ELI

## BONE LEVEL OCTA SYSTEM IMPRESSION COPING PICK UP



ø5.5

**HEX** S-ICOAP  
**NON-HEX** S-ICOAPN

## BONE LEVEL OCTA SYSTEM IMPRESSION COPING TRANSFER



ø5.0

**HEX** S-ICOAT  
**NON-HEX** S-ICOATN

**System** » Octa System  
**Material** » Ti Grade 5 ELI

## BONE LEVEL OCTA ABUTMENT LAB ANALOG



ø4.8

S-OALA

**System** » Octa System  
**Material** » Ti Grade 5 ELI

# BONE LEVEL O-RING SYSTEM

## BONE LEVEL O-RING ABUTMENTS



ø3.42



ø4.5



ø2.9\*



ø3.7\*

	ø3.42	ø4.5	ø2.9*	ø3.7*
<b>0.5 mm</b>	S-BORA05	-	S-BMORA05	-
<b>1 mm</b>	S-BORA10	-	-	-
<b>2 mm</b>	-	S-BORA20	-	S-BMORA20
<b>3 mm</b>	-	S-BORA30	-	-
<b>4 mm</b>	-	S-BORA40	-	S-BMORA40
<b>6 mm</b>	-	S-BORA60	-	-

\*only for ø3.3 implant

\*only for ø3.3 implant

**System** » Bone Level  
**Material** » Ti Grade 5 ELI  
**Torque Value** » 25 Ncm  
**Driver** » O-Ring Driver

## BONE LEVEL RETAINER



## BONE LEVEL O-RING LAB ANALOG



S-ORLA

**System** » Bone Level O-Ring System  
**Material** » Ti Grade 5 ELI

# LOCATOR SYSTEM

## BONE LEVEL MINI LOCATOR ABUTMENTS

**ø3.3**  
only for ø3.3 implant

<b>1.5 mm</b>	S-BMLA15
<b>2.5 mm</b>	S-BMLA25
<b>3.5 mm</b>	S-BMLA35
<b>4.5 mm</b>	S-BMLA45
<b>5.5 mm</b>	S-BMLA55
<b>6.5 mm</b>	S-BMLA65

## BONE LEVEL LOCATOR ABUTMENTS

**ø3.86**

<b>1.5 mm</b>	S-BLA15
<b>2.5 mm</b>	S-BLA25
<b>3.5 mm</b>	S-BLA35
<b>4.5 mm</b>	S-BLA45
<b>5.5 mm</b>	S-BLA55
<b>6.5 mm</b>	S-BLA65



**System** » Bone Level  
**Material** » Ti Grade 5 ELI  
**Torque Value** » 25 Ncm  
**Driver** » Locator Driver

## BONE LEVEL LOCATOR MATRIX



S-LMht



**VIOLET CAP**  
Rigid Retention 2.7 Kg



**YELLOW CAP**  
Extra Soft Retention 0.6 Kg



**CLEAR CAP**  
Standard Retention 1.8 Kg



**BLACK CAP**  
Processing



**PINK CAP**  
Soft Retention 1.2 Kg

## BONE LEVEL LOCATOR LAB ANALOG



S-BLLA Ø 4mm

# SUPERB OVERDENTURE SYSTEM

## SUPERB OVERDENTURE ABUTMENT

	ø2.9	ø3.5	ø3.42	ø4.5
<b>0.5 mm</b>	S-BMBA05	-	S-BEBA05	-
<b>1 mm</b>	-	-	S-BEBA10	-
<b>2 mm</b>	-	S-BMEBA20	-	S-BEBA20
<b>3 mm</b>	-	-	-	S-BEBA30
<b>4 mm</b>	-	S-BMEBA40	-	S-BEBA40
<b>5 mm</b>	-	-	-	S-BEBA50
<b>6 mm</b>	-	-	-	S-BEBA60



**System** » Bone Level  
**Material** » Ti Grade 5 ELI  
**Torque Value** » 25 Ncm  
**Driver** » Superb Overdenture Driver

## SUPERB OVERDENTURE CAP



S-EBOCap  
 ø 5.5mmx3.5mm

## SUPERB OVERDENTURE LAB ANALOG



S-EBLA  
 ø 4.1mmx16,6 mm

# MULTI UNIT SYSTEM

## MULTI UNIT ABUTMENTS



ø4.8 (17°)

**2.5 mm** S-BMUA48172  
**3.5 mm** S-BMUA48173



ø4.8 (35°)

S-BMUA48352  
 S-BMUA48353

**System** » Bone Level  
**Material** » Ti Grade 5 ELI  
**Torque Value** » : Abutment Screw  
 25 Ncm Oclusal Screw 20 Ncm  
**Driver** » 1,2 Hex Driver

## MULTI UNIT ABUTMENT HEALING CAP



ø4.8

S-BMUAHCap

**System** » Multi Unit System  
**Material** » Ti Grade 5 ELI

## MULTI UNIT IMPRESSION COPING TRANSFER



ø5.0

**NON-HEX** S-ICBMUATN

**System** » Multi Unit System  
**Material** » Ti Grade 5 ELI

## MULTI UNIT ABUTMENT PLASTIC CYLINDER



S-MUAPICY

## MULTI UNIT LAB ANALOG



S-MUALA  
 ø 4.8x14mm

**System** » Bone Level  
**Material** » Ti Grade 5 ELI



# BAR SYSTEM

**BAR BURNOUT** / S-IB - burnout

---



**BAR POSITIONING** / S-IBP - burnout

---



**BAR HOUSING** / S-IBH - burnout

---



**BAR CLIP YELLOW** / S-IBClip - y / Medium Retention 1.8 Kg

---



**BAR CLIP PINK** / S-IBClip - p / Soft Retention 1.2 Kg

---



**BAR CLIP RED** / S-IBClip - r / Extra Soft Retention 0.8 Kg

---



**TOOL FOR INSERTING CLIP**

---



















**KEY FOR PARALLELOMETER**

---



















# TISSUE LEVEL
















## COUPLE SYSTEM

FLOW	FIXTURE	HEALING	ABUTMENT	LAB ANALOG	ABUTMENT	SCREW
 <p>ø3.7    ø4.3-L06    ø4.3</p>  <p>ø4.8-L06    ø4.8    ø5.5</p>	 <p>Cover Screw</p>  <p>Healing Abutment</p>	<p>Pick-up</p>  <p>Hex    N-Hex</p> <p>Transfer</p>  <p>Hex    N-Hex</p>		<p>Couple</p>  <p>Hex    N-Hex</p> <p>Multi</p>  <p>Hex    N-Hex</p> <p>Angled</p>  <p>15°    25°</p> <p>Synocta</p>  <p>Hex    N-Hex</p>		
TOOL	 <p>Ratchet Driver</p>	 <p>Hex Screw Driver</p>	 <p>Hex Screw Driver</p>			 <p>Hex Screw Driver</p>












## SOLID SYSTEM

FLOW	FIXTURE	HEALING	ABUTMENT	IMPRESSION CAP	LAB ANALOG	PLASTIC COPING
 <p>ø3.7    ø4.3-L06    ø4.3</p>  <p>ø4.8-L06    ø4.8    ø5.5</p>	 <p>Closing Screw</p>  <p>Healing Abutment</p>	  <p>Cap</p>		<p>Shoulder Analog &amp; Pin</p>  <p>Lab Analog</p> 	<p>Crown</p>  <p>Bridge</p> 	
TOOL	 <p>Ratchet Driver</p>	 <p>Hex Screw Driver</p>	  <p>Solid Driver    Hex Screw Driver</p>			 <p>Hex Screw Driver</p>

## OCTA SYSTEM

FLOW	FIXTURE	HEALING	OCTA ABUTMENT	IMPRESSION COPING	LAB ANALOG	LAB CYLINDER	SCREW
 <p>ø3.7   ø4.3-L06   ø4.3</p>  <p>ø4.8-L06   ø4.8   ø5.5</p>	 <p>Cover Screw</p>  <p>Healing Abutment</p>	 <p>Octa Abutment</p>	 <p>Pick-up</p>  <p>Transfer</p>	 <p>Lab Analog</p>	 <p>Octa</p>  <p>N-Octa</p>	 <p>Screw</p>	
TOOL	 <p>Ratchet Driver</p>	 <p>Hex Screw Driver</p>	 <p>Octa Driver</p>				 <p>Hex Screw Driver</p>

## O-RING SYSTEM

FLOW	FIXTURE	HEALING	O-RING ABUTMENT	O-RING LAB ANALOG	RETAINER
 <p>ø3.7   ø4.3-L06   ø4.3</p>  <p>ø4.8-L06   ø4.8   ø5.5</p>	 <p>Closing Screw</p>  <p>Healing Abutment</p>	 <p>O-Ring Abutment</p>	 <p>O-Ring Lab Analog</p>	 <p>O-Ring</p>  <p>O-Ring Retainer</p>	
TOOL	 <p>Ratchet Driver</p>	 <p>Hex Screw Driver</p>	 <p>O-Ring Driver</p>		

# TISSUE LEVEL HEALING ABUTMENTS

## TISSUE LEVEL HEALING ABUTMENTS



ø3.7

ø6.4

<b>3 mm</b>	S-THA30	S-THAW30
<b>4 mm</b>	S-THA40	S-THAW40
<b>5.5 mm</b>	S-THA55	S-THAW55

**System** » Tissue Level  
**Material** » Ti Grade 4  
**Torque Value** » 20 Ncm  
**Driver** » Hex Screw Driver

# TISSUE LEVEL COUPLE SYSTEM

## TISSUE LEVEL COUPLE ABUTMENTS

## TISSUE LEVEL COUPLE ABUTMENT NON-OCTA



ø5.2

ø6.4

<b>1 mm</b>	S-TCA10	S-TCAW10
<b>2 mm</b>	S-TCA20	S-TCAW20
<b>3 mm</b>	S-TCA30	S-TCAW30
<b>4 mm</b>	S-TCA40	S-TCAW40

ø5.2

ø6.4

<b>1 mm</b>	S-TCAN10	S-TCAWN10
<b>2 mm</b>	S-TCAN20	S-TCAWN20
<b>3 mm</b>	S-TCAN30	S-TCAWN30
<b>4 mm</b>	S-TCAN40	S-TCAWN40

**System** » Tissue Level // **Material** » Ti Grade 4 // **Torque Value** » 25 Ncm // **Driver** » Hex Screw Driver

## TISSUE LEVEL ANGLED ABUTMENT



ø3.7

<b>15°</b>	S-TAA15
<b>25°</b>	S-TAA25

**System** » Tissue Level  
**Material** » Ti Grade 4  
**Torque Value** » 25 Ncm  
**Driver** » Hex Screw Driver

## TISSUE LEVEL MULTI ABUTMENT OCTA / NON-OCTA



ø5.2

<b>OCTA</b>	S-TMA
<b>N-OCTA</b>	S-TMAN

**System** » Tissue Level  
**Material** » Ti Grade 4  
**Torque Value** » 25 Ncm  
**Driver** » Hex Screw Driver

## TISSUE LEVEL SYNOCTA ABUTMENT



ø3.5

6 mm

S-TSOA

## TISSUE LEVEL SYNOCTA ABUTMENT NON-OCTA



ø3.5

6 mm

S-TSOAN

System » Tissue Level // Material » Ti Grade 4 // Torque Value » 25 Ncm // Driver » Hex Screw Driver

## TISSUE LEVEL IMPRESSION COPING PICK UP OCTA / NON-OCTA



ø5.5

OCTA

S-TICP

N-OCTA

S-TICPN



ø6.4

S-TICPW

S-TICPWN

System » Tissue Level  
Material » Ti Grade 5 ELI

## TISSUE LEVEL IMPRESSION COPING TRANSFER OCTA / NON-OCTA



ø5.5

OCTA

S-TICT

N-OCTA

S-TICTN

System » Tissue Level  
Material » Ti Grade 5 ELI

## TISSUE LEVEL FIXTURE LAB ANALOG



ø4.8

S-TFLA



ø6.0

S-TFLLA

System » Tissue Level  
Material » Ti Grade 5 ELI

# TISSUE LEVEL SOLID SYSTEM

## TISSUE LEVEL SOLID ABUTMENTS



ø3.5      ø4.5

<b>4 mm</b>	S-TSA40	S-TSAL40
<b>5 mm</b>	-	S-TSAL50
<b>5.5 mm</b>	S-TSA55	-
<b>7 mm</b>	S-TSA70	S-TSAL70

**System** » Tissue Level  
**Material** » Ti Grade 4  
**Torque Value** » 25 Ncm  
**Driver** » Solid Driver

## TISSUE LEVEL SOLID ABUTMENT POSITIONING CYLINDER



S-TSAPoCy

**System** » Tissue Level Solid System  
**Material** » Plastic

## TISSUE LEVEL SOLID PLASTIC COPING



ø5.3

<b>CROWN</b>	S-TSAPICc
<b>BRIDGE</b>	S-TSAPICb

**System** » Tissue Level Solid System  
**Material** » POM

## TISSUE LEVEL SOLID ABUTMENT IMPRESSION CAP



S-TSAICap

**System** » Tissue Level Solid System  
**Material** » Plastic

## TISSUE LEVEL SOLID ABUTMENT LAB ANALOG



ø4.8



ø6.0

<b>4 mm</b>	S-TSALA40	S-TSALLA40
<b>5.5 mm</b>	S-TSALA55	S-TSALLA55
<b>7 mm</b>	S-TSALA70	S-TSALLA70

**System** » Tissue Level Solid System  
**Material** » Ti Grade 5 ELI

## TISSUE LEVEL OCTA SYSTEM

### TISSUE LEVEL OCTA ABUTMENT



ø3.5

S-TOA

**System** » Tissue Level  
**Material** » Ti Grade 4  
**Torque Value** » : Abutment Screw  
 25 Ncm Oclusal Screw 20 Ncm  
**Driver** » Octa Driver

### TISSUE LEVEL OCTA ABUTMENT HEALING CAP



ø4.8

<b>4 mm</b>	S-OAHCap
<b>6 mm</b>	S-OAHCapL

**System** » Octa System  
**Material** » Ti Grade 5 ELI

### TISSUE LEVEL OCTA ABUTMENT PLASTIC CYLINDER OCTA / NON-OCTA



<b>OCTA</b>	S-OAPICy
<b>N-OCTA</b>	S-OAPICyN

**System** » Octa System  
**Material** » POM

## TISSUE LEVEL OCTA SYSTEM IMPRESSION COPING PICK UP OCTA / NON-OCTA



ø5.5

**OCTA** S-ICOAP  
**N-OCTA** S-ICOAPN

**System** » Octa System  
**Material** » Ti Grade 5 ELI

## TISSUE LEVEL OCTA SYSTEM IMPRESSION COPING TRANSFER OCTA / NON-OCTA



ø5.0

**OCTA** S-ICOAT  
**N-OCTA** S-ICOATN

**System** » Octa System  
**Material** » Ti Grade 5 ELI

## TISSUE LEVEL OCTA LAB ANALOG



ø4.8

S-OALA

**System** » Octa System  
**Material** » Ti Grade 5 ELI

# TISSUE LEVEL O-RING SYSTEM

## TISSUE LEVEL RETAINER



S-ORR



S-ORRc



TITANIUM CAP



GREEN RETAINER  
Retention 350 Gr



PINK RETAINER  
Retention 900 Gr



YELLOW RETAINER  
Retention 500 Gr

S-ORRo



## TISSUE LEVEL O-RING ABUTMENTS



Ø3.5

<b>0.5 mm</b>	S-TORA05
<b>2 mm</b>	S-TORA20
<b>4 mm</b>	S-TORA40

**System** » Tissue Level  
**Material** » Ti Grade 4  
**Torque Value** » 25 Ncm  
**Driver** » O-Ring Driver

## TISSUE LEVEL O-RING LAB ANALOG



Ø 4.1mm

S-ORLA

**System** » O-Ring System  
**Material** » Ti Grade 5 ELI

# TISSUE LEVEL LOCATOR SYSTEM

## TISSUE LEVEL LOCATOR ABUTMENTS

Ø4.9

<b>1.5 mm</b>	S-TLA15
<b>2.5 mm</b>	S-TLA25
<b>3.5 mm</b>	S-TLA35
<b>4.5 mm</b>	S-TLA45
<b>5.5 mm</b>	S-TLA55
<b>6.5 mm</b>	S-TLA65



**System** » Tissue Level  
**Material** » Ti Grade 5 ELI  
**Torque Value** » 25 Ncm  
**Driver** » Locator Driver

## TISSUE LEVEL LOCATOR MATRIX



S-LMht



**VIOLET CAP**  
Rigid Retention 2.7 Kg



**YELLOW CAP**  
Extra Soft Retention 0.6 Kg



**CLEAR CAP**  
Standard Retention 1.8 Kg



**BLACK CAP**  
Processing



**PINK CAP**  
Soft Retention 1.2 Kg

## TISSUE LEVEL LOCATOR LAB ANALOG



Ø 4.90mm

S-TLLA

# TISSUE LEVEL SUPERB OVERDENTURE SYSTEM

## TISSUE LEVEL SUPERB OVERDENTURE ABUTMENT

---

	<u>ø3.5</u>
<b>0.5 mm</b>	S-TEBA05
<b>2 mm</b>	S-TEBA20
<b>4 mm</b>	S-TEBA40

**System** » Bone Level  
**Material** » Ti Grade 5 ELI  
**Torque Value** » 25 Ncm  
**Driver** » Superb Overdenture Driver

## TISSUE LEVEL SUPERB OVERDENTURE CAP

---



## TISSUE LEVEL SUPERB OVERDENTURE LAB ANALOG

---



# ZYGOMA DENTAL IMPLANT SYSTEM

## ABUTMENTS

### ABUTMENTS



**3 mm** S-ZCA4330  
**5 mm** S-ZCA4350

### ANGLED ABUTMENTS



**3 mm** S-ZAA43203  
**5 mm** S-ZAA43205

**System** » Zygoma  
**Material** » Ti Grade 5 ELI  
**Torque Value** » 25 Ncm  
**Driver** » Hex Screw Driver

### MULTI UNIT



**17°** S-ZMUA43172

**System** » Zygoma  
**Material** » Ti Grade 5 ELI  
**Driver** » Hex Screw Driver

## LAB ANALOG



S-ZLA

**System** » Zygoma // **Material** » Ti Grade 5 ELI

## CAPS

### HEALING CAP



**3 mm** S-ZHA4330  
**5 mm** S-ZHA4350

**System** » Zygoma // **Torque Value** » 25 Ncm

## COPING

### IMPRESSION COPING PICK UP



S-ZICP

**System** » Zygoma  
**Material** » Ti Grade 5 ELI

## PLASTIC CYLINDER



S-ZPICYN

**System** » Zygoma  
**Material** » POM

# DRIVERS

## Implant Drivers (For ratchet and contra angle)

Bone Level Mini Fixture Driver Machine Short	S-BMFDMS
Bone Level Mini Fixture Driver Machine Long	S-BMFDML
Bone Level Mini Fixture Driver Ratchet Short	S-BMFDRS
Bone Level Mini Fixture Driver Ratchet Long	S-BMFDRL
Bone Level Fixture Driver Machine Short	S-BFDMS
Bone Level Fixture Driver Machine Long	S-BFDML
Bone Level Fixture Driver Ratchet Short	S-BFDRS
Bone Level Fixture Driver Ratchet Long	S-BFDRL
Tissue Level Fixture Driver Ratchet Short	S-TFDRS
Tissue Level Fixture Driver Ratchet Long	S-TFDRL
Tissue Level Fixture Driver Machine Short	S-TFDMS
Tissue Level Fixture Driver Machine Long	S-TFDML
Hex driver (extra short, short, long, extra long)	S-HDRS / S-HDRL
Hex Driver Contra Angle (extra short, short, long, extra long)	S-HDRS / S-HDRL
Superb Overdenture Driver	S-EBRD
O-ring Driver	S-ORADS/L
Octa Driver	S-OADS/L
Solid Driver	S-SADL
ORRO Insert Adaptor	S-OIA
Bar Matrix Insert Adaptor	S-BMIA

# SCREWS

## CLOSING SCREW

### BONE LEVEL

Bone Level Mini Closing Screw	ø 2.8	4.5 mm	S-BMCS
Bone Level Closing Screw	ø 3.36	6.5 mm	S-BCS

### TISSUE LEVEL

Tissue Level Closing-Screw Small	ø 3.5	5.35 mm	S-TCSS
Tissue Level Closing-Screw Large	ø 5.2	6.4mm	S-TCSL
Tissue Level Closing-Screw Wide	ø 6.4	6.4 mm	S-TCSW

## ABUTMENT SCREWS

### BONE LEVEL

Bone Level Mini Abutment Screw	ø 2.1	7.64 mm	S-BMAS
Bone Level Abutment Screw	ø 2.3	9.32 mm	S-BAS
Bone Level Multi-Unit Abutment Screw	ø 1.96	9 mm	S-BMUAS

### TISSUE LEVEL

Tissue Level Abutment Screw	ø 2.5	7.1	S-TAS
-----------------------------	-------	-----	-------

## COMMON SCREWS

Octa Abutment Screw	ø 2.5	4.5 mm	S-OAS
Multi-Unit Abutment Healing Screw	ø 2.2	4 mm	S-MUAHS

# INDEX CODE

PRODUCT	PLATFORM	HEIGHT	CODE
<b>BONE LEVEL FIXTURES</b>			
<b>BONE LEVEL MINI FIXTURE</b>			
Bone Level Mini Fixture Straight SRA	ø 3.3	8 mm	S-BMFSR3308
Bone Level Mini Fixture Straight SRA	ø 3.3	10 mm	S-BMFSR3310
Bone Level Mini Fixture Straight SRA	ø 3.3	12 mm	S-BMFSR3312
Bone Level Mini Fixture Straight SRA	ø 3.3	14 mm	S-BMFSR3314
<b>BONE LEVEL SHORT FIXTURE</b>			
Bone Level Short Fixture Straight SRA	ø 4.3	5 mm	S-BMFSR4305
Bone Level Short Fixture Straight SRA	ø 4.3	6 mm	S-BMFSR4306
Bone Level Short Fixture Straight SRA	ø 4.8	5 mm	S-BMFSR4805
Bone Level Short Fixture Straight SRA	ø 4.8	6 mm	S-BMFSR4806
Bone Level Short Fixture Straight SRA	ø 5.5	5 mm	S-BMFSR5505
Bone Level Short Fixture Straight SRA	ø 5.5	6 mm	S-BMFSR5506
<b>BONE LEVEL FIXTURE</b>			
Bone Level Wide Fixture Hybrid SRA	ø 3.7	8 mm	S-BFHR3708
Bone Level Wide Fixture Hybrid SRA	ø 3.7	10 mm	S-BFHR3710
Bone Level Wide Fixture Hybrid SRA	ø 3.7	12 mm	S-BFHR3712
Bone Level Wide Fixture Hybrid SRA	ø 3.7	14 mm	S-BFHR3714
Bone Level Wide Fixture Hybrid SRA	ø 4.3	8 mm	S-BFHR4308
Bone Level Wide Fixture Hybrid SRA	ø 4.3	10 mm	S-BFHR4310
Bone Level Wide Fixture Hybrid SRA	ø 4.3	12 mm	S-BFHR4312
Bone Level Wide Fixture Hybrid SRA	ø 4.3	14 mm	S-BFHR4314
Bone Level Wide Fixture Hybrid SRA	ø 4.8	8 mm	S-BFHR4808
Bone Level Wide Fixture Hybrid SRA	ø 4.8	10 mm	S-BFHR4810
Bone Level Wide Fixture Hybrid SRA	ø 4.8	12 mm	S-BFHR4812
Bone Level Wide Fixture Hybrid SRA	ø 4.8	14 mm	S-BFHR4814
<b>BONE LEVEL WIDE FIXTURE</b>			
Bone Level Wide Fixture Hybrid SRA	ø 5.5	8 mm	S-BFHR5508
Bone Level Wide Fixture Hybrid SRA	ø 5.5	10 mm	S-BFHR5510
Bone Level Wide Fixture Hybrid SRA	ø 5.5	12 mm	S-BFHR5512
<b>HEALING ABUTMENTS</b>			
<b>HEALING ABUTMENTS</b>			
Healing Abutment	ø 3.5	2.5 mm	S-BMHA3725
Healing Abutment	ø 3.5	3.5 mm	S-BMHA3735
Healing Abutment	ø 3.5	6 mm	S-BMHA3750
Healing Abutment	ø 3.5	7 mm	S-BMHA3770
Healing Abutment	ø 4.5	1 mm	S-BHA4510
Healing Abutment	ø 4.5	2 mm	S-BHA4520

PRODUCT	PLATFORM	HEIGHT	CODE
Healing Abutment	∅ 4.5	3 mm	S-BHA4530
Healing Abutment	∅ 4.5	4 mm	S-BHA4540
Healing Abutment	∅ 4.5	5 mm	S-BHA4550
Healing Abutment	∅ 5.5	1 mm	S-BHA5510
Healing Abutment	∅ 5.5	2 mm	S-BHA5520
Healing Abutment	∅ 5.5	3 mm	S-BHA5530
Healing Abutment	∅ 5.5	4 mm	S-BHA5540
Healing Abutment	∅ 5.5	5 mm	S-BHA5550
Healing Abutment	∅ 6.5	1 mm	S-BHA6510
Healing Abutment	∅ 6.5	2 mm	S-BHA6520
Healing Abutment	∅ 6.5	3 mm	S-BHA6530
Healing Abutment	∅ 6.5	4 mm	S-BHA6540
Healing Abutment	∅ 6.5	5 mm	S-BHA6550
Mini Healing Abutment for Short Implant	∅ 5.0	2 mm	S-BMHA5020
Mini Healing Abutment for Short Implant	∅ 5.0	3 mm	S-BMHA5030
Mini Healing Abutment for Short Implant	∅ 6.0	2 mm	S-BMHA6020
Mini Healing Abutment for Short Implant	∅ 6.0	3 mm	S-BMHA6030
Mini Healing Abutment for Short Implant	∅ 8.0	2 mm	S-BMHA8020
Mini Healing Abutment for Short Implant	∅ 8.0	3 mm	S-BMHA8030

## COUPLE SYSTEM

### COUPLE ABUTMENTS

Mini Couple Abutment H7	∅ 3.5	1 mm	S-BMCA37107
Mini Couple Abutment H7	∅ 3.5	2 mm	S-BMCA37207
Mini Couple Abutment H7	∅ 3.5	3 mm	S-BMCA37307
Mini Couple Abutment H7	∅ 3.5	4 mm	S-BMCA37407
Mini Couple Abutment H7	∅ 3.5	5 mm	S-BMCA37507
Couple Abutment H7	∅ 4.5	1 mm	S-BCA45107
Couple Abutment H7	∅ 4.5	2 mm	S-BCA45207
Couple Abutment H7	∅ 4.5	3 mm	S-BCA45307
Couple Abutment H7	∅ 4.5	4 mm	S-BCA45407
Couple Abutment H7	∅ 4.5	5 mm	S-BCA45507
Couple Abutment H7	∅ 5.5	1 mm	S-BCA55107
Couple Abutment H7	∅ 5.5	2 mm	S-BCA55207
Couple Abutment H7	∅ 5.5	3 mm	S-BCA55307
Couple Abutment H7	∅ 5.5	4 mm	S-BCA55407
Couple Abutment H7	∅ 5.5	5 mm	S-BCA55507
Couple Abutment H7	∅ 6.5	1 mm	S-BCA65107
Couple Abutment H7	∅ 6.5	2 mm	S-BCA65207
Couple Abutment H7	∅ 6.5	3 mm	S-BCA65307
Couple Abutment H7	∅ 6.5	4 mm	S-BCA65407
Couple Abutment H7	∅ 6.5	5 mm	S-BCA65507

PRODUCT	PLATFORM	HEIGHT	CODE
Mini Couple Abutment for Short Implant	ø 5.0	2 mm	S-BMCA5020
Mini Couple Abutment for Short Implant	ø 5.0	3 mm	S-BMCA5030
Mini Couple Abutment for Short Implant	ø 6.0	2 mm	S-BMCA6020
Mini Couple Abutment for Short Implant	ø 6.0	3 mm	S-BMCA6030
Mini Couple Abutment for Short Implant	ø 8.0	2 mm	S-BMCA8020
Mini Couple Abutment for Short Implant	ø 8.0	3 mm	S-BMCA8030
<b>COUPLE ABUTMENTS NON-HEX</b>			
Mini Couple Abutment H7 Non-Hex	ø 3.5	1 mm	S-BMCAN3710
Mini Couple Abutment H7 Non-Hex	ø 3.5	2 mm	S-BMCAN3720
Mini Couple Abutment H7 Non-Hex	ø 3.5	3 mm	S-BMCAN3730
Mini Couple Abutment H7 Non-Hex	ø 3.5	4 mm	S-BMCAN3740
Mini Couple Abutment H7 Non-Hex	ø 3.5	5 mm	S-BMCAN3750
Couple Abutment H7 Non-Hex	ø 4.5	1 mm	S-BCAN4510
Couple Abutment H7 Non-Hex	ø 4.5	2 mm	S-BCAN4520
Couple Abutment H7 Non-Hex	ø 4.5	3 mm	S-BCAN4530
Couple Abutment H7 Non-Hex	ø 4.5	4 mm	S-BCAN4540
Couple Abutment H7 Non-Hex	ø 4.5	5 mm	S-BCAN4550
Couple Abutment H7 Non-Hex	ø 5.5	1 mm	S-BCAN5510
Couple Abutment H7 Non-Hex	ø 5.5	2 mm	S-BCAN5520
Couple Abutment H7 Non-Hex	ø 5.5	3 mm	S-BCAN5530
Couple Abutment H7 Non-Hex	ø 5.5	4 mm	S-BCAN5540
Couple Abutment H7 Non-Hex	ø 5.5	5 mm	S-BCAN5550
Couple Abutment H7 Non-Hex	ø 6.5	1 mm	S-BCAN6510
Couple Abutment H7 Non-Hex	ø 6.5	2 mm	S-BCAN6520
Couple Abutment H7 Non-Hex	ø 6.5	3 mm	S-BCAN6530
Couple Abutment H7 Non-Hex	ø 6.5	4 mm	S-BCAN6540
Couple Abutment H7 Non-Hex	ø 6.5	5 mm	S-BCAN6550
<b>ANGLED ABUTMENTS</b>			
Mini Angled Abutment	ø 3.5	1 mm 15°	S-BMAA37151
Mini Angled Abutment	ø 3.5	2 mm 15°	S-BMAA37152
Mini Angled Abutment	ø 3.5	4 mm 15°	S-BMAA37154
Angled Abutment	ø 4.5	2 mm 15°	S-BAA45152
Angled Abutment	ø 4.5	2 mm 25°	S-BAA45252
Angled Abutment	ø 4.5	4 mm 15°	S-BAA45154
Angled Abutment	ø 4.5	4 mm 25°	S-BAA45254
Angled Abutment	ø 5.5	2 mm 15°	S-BAA55152
Angled Abutment	ø 5.5	2 mm 25°	S-BAA55252
Angled Abutment	ø 5.5	4 mm 15°	S-BAA55154
Angled Abutment	ø 5.5	4 mm 25°	S-BAA55254
Angled Abutment	ø 6.5	2 mm 15°	S-BAA65152
Angled Abutment	ø 6.5	2 mm 25°	S-BAA65252
Angled Abutment	ø 6.5	4 mm 15°	S-BAA65154
Angled Abutment	ø 6.5	4 mm 25°	S-BAA65254



<b>PRODUCT</b>	<b>PLATFORM</b>	<b>HEIGHT</b>	<b>CODE</b>
<b>ANGLED ABUTMENTS NON-HEX</b>			
Mini Angled Abutment Non-Hex	ø 3.5	1 mm 15°	S-BMAAN37151
Mini Angled Abutment Non-Hex	ø 3.5	2 mm 15°	S-BMAAN37152
Mini Angled Abutment Non-Hex	ø 3.5	4 mm 15°	S-BMAAN37154
Angled Abutment Non-Hex	ø 4.5	2 mm 15°	S-BAAN45152
Angled Abutment Non-Hex	ø 4.5	2 mm 25°	S-BAAN45252
Angled Abutment Non-Hex	ø 4.5	4 mm 15°	S-BAAN45154
Angled Abutment Non-Hex	ø 4.5	4 mm 25°	S-BAAN45254
Angled Abutment Non-Hex	ø 5.5	2 mm 15°	S-BAAN55152
Angled Abutment Non-Hex	ø 5.5	2 mm 25°	S-BAAN55252
Angled Abutment Non-Hex	ø 5.5	4 mm 15°	S-BAAN55154
Angled Abutment Non-Hex	ø 5.5	4 mm 25°	S-BAAN55254
Angled Abutment Non-Hex	ø 6.5	2 mm 15°	S-BAAN65152
Angled Abutment Non-Hex	ø 6.5	2 mm 25°	S-BAAN65252
Angled Abutment Non-Hex	ø 6.5	4 mm 15°	S-BAAN65154
Angled Abutment Non-Hex	ø 6.5	4 mm 25°	S-BAAN65254
<b>MULTI ABUTMENTS</b>			
Mini Multi Abutment	ø 3.5		S-BMMA37
Multi Abutment	ø 4.5		S-BMA45
Multi Abutment	ø 5.5		S-BMA55
<b>MULTI ABUTMENTS NON-HEX</b>			
Mini Multi Abutment Non-Hex	ø 3.5		S-BMMAN37
Multi Abutment Non-Hex	ø 4.5		S-BMAN45
Multi Abutment Non-Hex	ø 5.5		S-BMAN55
<b>TEMPORARY ABUTMENTS</b>			
Mini Temporary Abutment	ø 3.7		S-BMTA37
Temporary Abutment	ø 4.5		S-BTA45
Temporary Abutment	ø 5.5		S-BTA55
<b>IMPRESSION COPING PICK UP</b>			
Mini Impression Coping Pick up	ø 3.5		S-BMICP37L
Mini Impression Coping Pick up	ø 3.5		S-BMICP37S
Impression Coping Pick up	ø 4.5		S-BICP45L
Impression Coping Pick up	ø 4.5		S-BICP45S
Impression Coping Pick up	ø 5.5		S-BICP55L
Impression Coping Pick up	ø 5.5		S-BICP55S
Impression Coping Pick up	ø 6.5		S-BICP65L
Impression Coping Pick up	ø 6.5		S-BICP65S
<b>IMPRESSION COPING PICK UP NON-HEX</b>			
Mini Impression Coping Pick up Non-Hex	ø 3.5		S-BMICPN37L
Mini Impression Coping Pick up Non-Hex	ø 3.5		S-BMICPN37S

<b>PRODUCT</b>	<b>PLATFORM</b>	<b>HEIGHT</b>	<b>CODE</b>
Impression Coping Pick up Non-Hex	ø 4.5		S-BICPN45L
Impression Coping Pick up Non-Hex	ø 4.5		S-BICPN45S
Impression Coping Pick up Non-Hex	ø 5.5		S-BICPN55L
Impression Coping Pick up Non-Hex	ø 5.5		S-BICPN55S
Impression Coping Pick up Non-Hex	ø 6.5		S-BICPN65L
Impression Coping Pick up Non-Hex	ø 6.5		S-BICPN65S
<b>IMPRESSION COPING TRANSFER</b>			
Mini Impression Coping Transfer	ø 3.5		S-BMICT37L
Mini Impression Coping Transfer	ø 3.5		S-BMICT37S
Impression Coping Transfer	ø 4.5		S-BICT45L
Impression Coping Transfer	ø 4.5		S-BICT45S
Impression Coping Transfer	ø 5.5		S-BICT55L
Impression Coping Transfer	ø 5.5		S-BICT55S
Impression Coping Transfer	ø 6.5		S-BICT65L
Impression Coping Transfer	ø 6.5		S-BICT65S
<b>IMPRESSION COPING TRANSFER NON-HEX</b>			
Mini Impression Coping Transfer Non-Hex	ø 3.5		S-BMICTN37L
Mini Impression Coping Transfer Non-Hex	ø 3.5		S-BMICTN37S
Impression Coping Transfer Non-Hex	ø 4.5		S-BICTN45L
Impression Coping Transfer Non-Hex	ø 4.5		S-BICTN45S
Impression Coping Transfer Non-Hex	ø 5.5		S-BICTN55L
Impression Coping Transfer Non-Hex	ø 5.5		S-BICTN55S
Impression Coping Transfer Non-Hex	ø 6.5		S-BICTN65L
Impression Coping Transfer Non-Hex	ø 6.5		S-BICTN65S
<b>FIXTURE LAB ANALOG</b>			
Mini Fixture Lab Analog	ø 3.3	12 mm	S-BMFLA
Fixture Lab Analog	ø 4.3	12 mm	S-BFLA
<b>SOLID SYSTEM</b>			
<b>SOLID ABUTMENTS</b>			
Mini Solid Abutment	ø 3.5	1 mm	S-BMSA3710
Mini Solid Abutment	ø 3.5	2 mm	S-BMSA3720
Mini Solid Abutment	ø 3.5	3 mm	S-BMSA3730
Mini Solid Abutment	ø 3.5	4 mm	S-BMSA3740
Mini Solid Abutment	ø 3.5	5 mm	S-BMSA3750

<b>PRODUCT</b>	<b>PLATFORM</b>	<b>HEIGHT</b>	<b>CODE</b>
Solid Abutment	∅ 4.5	1 mm	S-BSA4510
Solid Abutment	∅ 4.5	2 mm	S-BSA4520
Solid Abutment	∅ 4.5	3 mm	S-BSA4530
Solid Abutment	∅ 4.5	4 mm	S-BSA4540
Solid Abutment	∅ 4.5	5 mm	S-BSA4550
Solid Abutment	∅ 5.5	1 mm	S-BSA5510
Solid Abutment	∅ 5.5	2 mm	S-BSA5520
Solid Abutment	∅ 5.5	3 mm	S-BSA5530
Solid Abutment	∅ 5.5	4 mm	S-BSA5540
Solid Abutment	∅ 5.5	5 mm	S-BSA5550
Solid Abutment	∅ 6.5	1 mm	S-BSA6510
Solid Abutment	∅ 6.5	2 mm	S-BSA6520
Solid Abutment	∅ 6.5	3 mm	S-BSA6530
Solid Abutment	∅ 6.5	4 mm	S-BSA6540
Solid Abutment	∅ 6.5	5 mm	S-BSA6550

### **SOLID PLASTIC BURN-OUT COPING**

Mini Solid Plastic Coping Bridge	∅ 3.5		S-BMSAPIC37b
Solid Plastic Coping Bridge	∅ 4.5		S-BSAPIC45b
Solid Plastic Coping Bridge	∅ 5.5		S-BSAPIC55b
Solid Plastic Coping Bridge	∅ 6.5		S-BSAPIC65b
Mini Solid Plastic Coping Crown	∅ 3.5		S-BMSAPIC37c
Solid Plastic Coping Crown	∅ 4.5		S-BSAPIC45c
Solid Plastic Coping Crown	∅ 5.5		S-BSAPIC55c
Solid Plastic Coping Crown	∅ 6.5		S-BSAPIC65c

### **SOLID ABUTMENT IMPRESSION CAPS**

Mini Solid Abutment Impression Cap	∅ 3.5		S-BMSAICap37
Solid Abutment Impression Cap	∅ 4.5		S-BSAICap45
Solid Abutment Impression Cap	∅ 5.5		S-BSAICap55
Solid Abutment Impression Cap	∅ 6.5		S-BSAICap65

### **SOLID LAB ANALOG**

Mini Solid Lab Analog	∅ 3.5		S-BMSALA37
Solid Lab Analog	∅ 4.5		S-BSALA45
Solid Lab Analog	∅ 5.5		S-BSALA55
Solid Lab Analog	∅ 6.5		S-BSALA65

### **OCTA SYSTEM**

#### **OCTA ABUTMENTS**

Mini Octa Abutment	∅ 3.5	1 mm	S-BMOA10
Mini Octa Abutment	∅ 3.5	2 mm	S-BMOA20
Mini Octa Abutment	∅ 3.5	3 mm	S-BMOA30
Mini Octa Abutment	∅ 3.5	4 mm	S-BMOA40
Mini Octa Abutment	∅ 3.5	5 mm	S-BMOA50

<b>PRODUCT</b>	<b>PLATFORM</b>	<b>HEIGHT</b>	<b>CODE</b>
Octa Abutment	ø 4.8	1 mm	S-BOA10
Octa Abutment	ø 4.8	2 mm	S-BOA20
Octa Abutment	ø 4.8	3 mm	S-BOA30
Octa Abutment	ø 4.8	4 mm	S-BOA40
Octa Abutment	ø 4.8	5 mm	S-BOA50
<b>OCTA ABUTMENT HEALING CAP</b>			
Octa Abutment Healing Cap	ø 4.9	4 mm	S-OAHCap
Octa Abutment Healing Cap Large	ø 4.9	6 mm	S-OAHCAP-L
<b>OCTA SYSTEM PLASTIC CYLINDER OCTA</b>			
Octa System Plastic Cylinder Octa	ø 5.2 mm		S-OAPICy
<b>OCTA SYSTEM PLASTIC CYLINDER NON-OCTA</b>			
Octa System Plastic Cylinder Non-Octa	ø 5.2 mm		S-OAPICyN
<b>OCTA SYSTEM TITANIUM CYLINDER</b>			
Octa System Plastic Cylinder	ø 5.2 mm		S-OATICy
<b>OCTA SYSTEM TITANIUM CYLINDER NON-HEX</b>			
Octa System Plastic Cylinder Non-Hex	ø 5.2		S-OATICyN
<b>OCTA SYSTEM IMPRESSION COPING PICK UP</b>			
Octa System Impression Coping Pick Up	ø 5.5	10 mm	S-ICOAP
<b>OCTA SYSTEM IMPRESSION COPING PICK UP NON-HEX</b>			
Octa System Impression Coping Pick Up Non-Hex	ø 5.5	10 mm	S-ICOAPN
<b>OCTA SYSTEM IMPRESSION COPING TRANSFER</b>			
Octa System Impression Coping Transfer	ø 5.0	7 mm	S-ICOAT
<b>OCTA SYSTEM IMPRESSION COPING TRANSFER NON-HEX</b>			
Octa System Impression Coping Transfer Non-Hex	ø 5.0	7 mm	S-ICOATN
<b>OCTA ABUTMENT LAB ANALOG</b>			
Octa Abutment Lab Analog	ø 4.8	14 mm	S-OALA
<b>O-RING SYSTEM</b>			
<b>O-RING ABUTMENTS</b>			
Mini O-Ring Abutment	ø 2.9	0.5 mm	S-BMORA05
Mini O-Ring Abutment	ø 3.5	2 mm	S-BMORA20
Mini O-Ring Abutment	ø 3.5	4 mm	S-BMORA40
O-Ring Abutment	ø 3.42	0.5 mm	S-BORA05
O-Ring Abutment	ø 3.42	1 mm	S-BORA10
O-Ring Abutment	ø 4.5	2 mm	S-BORA20
O-Ring Abutment	ø 4.5	3 mm	S-BORA30
O-Ring Abutment	ø 4.5	4 mm	S-BORA40
O-Ring Abutment	ø 4.5	6 mm	S-BORA60

<b>PRODUCT</b>	<b>PLATFORM</b>	<b>HEIGHT</b>	<b>CODE</b>
<b>RETAINERS</b>			
Retainer	ø 5.5	2 mm	S-ORR
Retainer	ø 5.5	3.5 mm	S-ORR-C
Retainer	ø 5.5	2 mm	S-ORR-O
<b>O-RING LAB ANALOG</b>			
O-Ring Lab Analog	ø 4.1		S-ORLA
<b>LOCATOR SYSTEM</b>			
<b>LOCATOR ABUTMENTS</b>			
Locator Mini Abutment	ø 3.86mm	1.5mm	S-BMLA15
Locator Mini Abutment	ø 3.86mm	2.5mm	S-BMLA25
Locator Mini Abutment	ø 3.86mm	3.5mm	S-BMLA35
Locator Mini Abutment	ø 3.86mm	4.5mm	S-BMLA45
Locator Mini Abutment	ø 3.86mm	5.5mm	S-BMLA55
Locator Mini Abutment	ø 3.86mm	6.5mm	S-BMLA65
Locator Abutment	ø 3.86mm	1.5mm	S-BLA15
Locator Abutment	ø 3.86mm	2.5mm	S-BLA25
Locator Abutment	ø 3.86mm	3.5mm	S-BLA35
Locator Abutment	ø 3.86mm	4.5mm	S-BLA45
Locator Abutment	ø 3.86mm	5.5mm	S-BLA55
Locator Abutment	ø 3.86mm	6.5mm	S-BLA65
<b>LOCATOR MATRIX</b>			
Locator Matrix (includes processing cap)	ø 5.5mm		S-LMHt
Locator Cap (includes 4 pcs caps)	ø 4.86 mm		S-LRI
<b>LOCATOR ANALOG</b>			
Locator Lab Analog	ø 4 mm		S-BLLA
<b>SUPERB OVERDENTURE SYSTEM</b>			
<b>SUPERB OVERDENTURE ABUTMENT</b>			
Superb Overdenture Mini Abutment	ø 2.9 mm	0.5 mm	S-BMEBA05
Superb Overdenture Mini Abutment	ø 3,5 mm	20 mm	S-BMEBA20
Superb Overdenture Mini Abutment	ø 3.5mm	40 mm	S-BMEBA40
Superb Overdenture Abutment	ø 3.42 mm	0.5 mm	S-BEBA05
Superb Overdenture Abutment	ø 3.42 mm	1 mm	S-BEBA10
Superb Overdenture Abutment	ø 4.50 mm	2 mm	S-BEBA20
Superb Overdenture Abutment	ø 4.50 mm	3 mm	S-BEBA30
Superb Overdenture Abutment	ø 4.50 mm	4 mm	S-BEBA40
Superb Overdenture Abutment	ø 4.50 mm	5 mm	S-BEBA50
Superb Overdenture Abutment	ø 4.50 mm	6 mm	S-BEBA60

<b>PRODUCT</b>	<b>PLATFORM</b>	<b>HEIGHT</b>	<b>CODE</b>
<b>SUPERB OVERDENTURE CAP</b>			
Superb Overdenture Cap	ø 5.5 mm	3.5 mm	S-EBOCAP
<b>SUPERB OVERDENTURE LAB ANALOG</b>			
Superb Overdenture Analog	ø 4.1mm	16.6 mm	S-EBLA
<b>MULTI-UNIT SYSTEM</b>			
<b>MULTI-UNIT ABUTMENT</b>			
Bone Level Multi Unit Abutment	ø 4.8x17	2.5 mm	S-BMUA48172
Bone Level Multi Unit Abutment	ø 4.8x17	3.5 mm	S-BMUA48173
Bone Level Multi Unit Abutment	ø 4.8x35	2.5 mm	S-BMUA48352
Bone Level Multi Unit Abutment	ø 4.8x35	3.5 mm	S-BMUA48353
<b>MULTI UNIT HEALING CAP</b>			
Multi Unit Healing Cap	ø 4.9	4.0 mm	S-BMUAHCAP
<b>MULTI UNIT IMPRESSION COPING PICK UP</b>			
Bone Level Multi Unit Impression Coping Pick Up	ø 4.8		S-BMUAICP
<b>MULTI UNIT ABUTMENT PLASTIC CYLINDER</b>			
Multi Unit Abutment Plastic Cylinder	ø 4.9		S-MUAPICY
<b>MULTI UNIT LAB ANALOG</b>			
Multi Unit Abutment Lab Analog	ø 4.8	14 mm	S-BMUALA
<b>BAR SYSTEM</b>			
Bar burnout			S-IB-BURNOUT
Bar Positioning			S-IBP-BURNOUT
Bar Housing			S-IBH-BURNOUT
Bar Clip Yellow			S-IBClip-y
Bar Clip Pink			S-IBClip-P
Bar Clip Red			S-IBClip-r
Tool for Inserting Clip			
Key for Parallelometer			
<b>TISSUE LEVEL FIXTURES</b>			
<b>TISSUE LEVEL FIXTURE HYBRID</b>			
Tissue Level Wide Fixture Hybrid SRA	ø 3.7	8 mm	S-TFHR3708
Tissue Level Wide Fixture Hybrid SRA	ø 3.7	10 mm	S-TFHR3710
Tissue Level Wide Fixture Hybrid SRA	ø 3.7	12 mm	S-TFHR3712
Tissue Level Wide Fixture Hybrid SRA	ø 3.7	14 mm	S-TFHR3714
Tissue Level Wide Fixture Hybrid SRA	ø 4.3	8 mm	S-TFHR4308
Tissue Level Wide Fixture Hybrid SRA	ø 4.3	10 mm	S-TFHR4310
Tissue Level Wide Fixture Hybrid SRA	ø 4.3	12 mm	S-TFHR4312
Tissue Level Wide Fixture Hybrid SRA	ø 4.3	14 mm	S-TFHR4314

<b>PRODUCT</b>	<b>PLATFORM</b>	<b>HEIGHT</b>	<b>CODE</b>
Tissue Level Wide Fixture Hybrid SRA	ø 4.8	8 mm	S-TFHR4808
Tissue Level Wide Fixture Hybrid SRA	ø 4.8	10 mm	S-TFHR4810
Tissue Level Wide Fixture Hybrid SRA	ø 4.8	12 mm	S-TFHR4812
Tissue Level Wide Fixture Hybrid SRA	ø 4.8	14 mm	S-TFHR4814
<b>TISSUE LEVEL SHORT FIXTURE STRAIGHT</b>			
Tissue Level Short Fixture Straight SRA	ø 4.3	4 mm	S-TFSR4304
Tissue Level Short Fixture Straight SRA	ø 4.3	5 mm	S-TFSR4305
Tissue Level Short Fixture Straight SRA	ø 4.3	6 mm	S-TFSR4306
Tissue Level Short Fixture Straight SRA	ø 4.8	4 mm	S-TFSR4804
Tissue Level Short Fixture Straight SRA	ø 4.8	5 mm	S-TFSR4805
Tissue Level Short Fixture Straight SRA	ø 4.8	6 mm	S-TFSR4806
Tissue Level Short Fixture Straight SRA	ø 5.5	4 mm	S-TWFSR5504
Tissue Level Short Fixture Straight SRA	ø 5.5	5 mm	S-TWFSR5505
Tissue Level Short Fixture Straight SRA	ø 5.5	6 mm	S-TWFSR5506
<b>TISSUE LEVEL WIDE FIXTURE HYBRID</b>			
Tissue Level Wide Fixture Hybrid SRA	ø 5.5	8 mm	S-BFHR5508
Tissue Level Wide Fixture Hybrid SRA	ø 5.5	10 mm	S-BFHR5510
Tissue Level Wide Fixture Hybrid SRA	ø 5.5	12 mm	S-BFHR5512
<b>TISSUE LEVEL HEALING ABUTMENTS</b>			
<b>TISSUE LEVEL HEALING ABUTMENTS</b>			
Tissue Level Healing Abutment	ø 3.7	3 mm	S-THA30
Tissue Level Healing Abutment	ø 3.7	4 mm	S-THA40
Tissue Level Healing Abutment	ø 3.7	5.5 mm	S-THA55
Tissue Level Healing Abutment	ø 6.4	3 mm	S-THAW30
Tissue Level Healing Abutment	ø 6.4	4 mm	S-THAW40
Tissue Level Healing Abutment	ø 6.4	5.5 mm	S-THAW55
<b>TISSUE LEVEL COUPLE SYSTEM</b>			
<b>TISSUE LEVEL COUPLE ABUTMENTS</b>			
Tissue Level Couple Abutment	ø 5.2	1 mm	S-TCA10
Tissue Level Couple Abutment	ø 5.2	2 mm	S-TCA20
Tissue Level Couple Abutment	ø 5.2	3 mm	S-TCA30
Tissue Level Couple Abutment	ø 5.2	4 mm	S-TCA40
Tissue Level Couple Abutment Wide	ø 6.4	1 mm	S-TCAW10
Tissue Level Couple Abutment Wide	ø 6.4	2 mm	S-TCAW20
Tissue Level Couple Abutment Wide	ø 6.4	3 mm	S-TCAW30
Tissue Level Couple Abutment Wide	ø 6.4	4 mm	S-TCAW40

<b>PRODUCT</b>	<b>PLATFORM</b>	<b>HEIGHT</b>	<b>CODE</b>
<b>TISSUE LEVEL COUPLE ABUTMENT NON-OCTA</b>			
Tissue Level Couple Abutment Non-Octa	ø 5.2	1 mm	S-TCAN10
Tissue Level Couple Abutment Non-Octa	ø 5.2	2 mm	S-TCAN20
Tissue Level Couple Abutment Non-Octa	ø 5.2	3 mm	S-TCAN30
Tissue Level Couple Abutment Non-Octa	ø 5.2	4 mm	S-TCAN40
<b>TISSUE LEVEL COUPLE ABUTMENT WIDE NON-OCTA</b>			
Tissue Level Couple Abutment Non-Octa	ø 6.4	1 mm	S-TCAWN10
Tissue Level Couple Abutment Non-Octa	ø 6.4	2 mm	S-TCAWN20
Tissue Level Couple Abutment Non-Octa	ø 6.4	3 mm	S-TCAWN30
Tissue Level Couple Abutment Non-Octa	ø 6.4	4 mm	S-TCAWN40
<b>TISSUE LEVEL ANGLED ABUTMENTS</b>			
Tissue Level Angled Abutment	ø 3.7	15°	S-TAA15
Tissue Level Angled Abutment	ø 3.7	25°	S-TAA25
<b>TISSUE LEVEL MULTI ABUTMENT OCTA</b>			
Tissue Level Multi Abutment Octa	ø 5.2	1.3 mm	S-TMA
<b>TISSUE LEVEL MULTI ABUTMENT NON-OCTA</b>			
Tissue Level Multi Abutment Non-Octa	ø 5.2	1.3 mm	S-TMAN
<b>TISSUE LEVEL SYNOCTA ABUTMENT OCTA</b>			
Tissue Level Synocta Abutment Octa	ø 3.5	6 mm	S-TSOA
<b>TISSUE LEVEL SYNOCTA ABUTMENT NON-OCTA</b>			
Tissue Level Synocta Abutment Non-Octa	ø 3.5	6 mm	S-TSOAN
<b>TISSUE LEVEL IMPRESSION COPING PICK UP OCTA</b>			
Tissue Level Impression Coping Pick Up Octa	ø 5.5	11.75 mm	S-TICP
Tissue Level Impression Coping Pick Up Wide Octa	ø 6.4	11.75 mm	S-TICPW
<b>TISSUE LEVEL IMPRESSION COPING PICK UP NON-OCTA</b>			
Tissue Level Impression Coping Pick Up Non-Octa	ø 5.5	11.75 mm	S-TICPN
Tissue Level Impression Coping Pick Up Wide Non-Octa	ø 6.4	11.75 mm	S-TICPWN
<b>TISSUE LEVEL IMPRESSION COPING TRANSFER</b>			
Tissue Level Impression Coping Transfer Octa	ø 5.5	9 mm	S-TICT
<b>TISSUE LEVEL IMPRESSION COPING TRANSFER NON-OCTA</b>			
Tissue Level Impression Coping Transfer Non-Octa	ø 5.5	9 mm	S-TICTN
<b>TISSUE LEVEL FIXTURE LAB ANALOG</b>			
Tissue Level Fixture Lab Analog	ø 4.8	14 mm	S-TFLA
Tissue Level Fixture Large Lab Analog	ø 6.0	14 mm	S-TFLLA



PRODUCT	PLATFORM	HEIGHT	CODE
<b>TISSUE LEVEL SOLID SYSTEM</b>			
<b>TISSUE LEVEL SOLID ABUTMENTS</b>			
Tissue Level Solid Abutment	ø 3.5	4 mm	S-TSA40
Tissue Level Solid Abutment	ø 3.5	5.5 mm	S-TSA55
Tissue Level Solid Abutment	ø 3.5	7 mm	S-TSA70
Tissue Level Solid Abutment Large	ø 4.5	4 mm	S-TSAL40
Tissue Level Solid Abutment Large	ø 4.5	5 mm	S-TSAL50
Tissue Level Solid Abutment Large	ø 4.5	7 mm	S-TSAL70
<b>TISSUE LEVEL SOLID ABUTMENT POSITIONING CYLINDER</b>			
Tissue Level Solid Abutment Positioning Cylinder	ø 10.2		S-TSAPoCy
<b>TISSUE LEVEL SOLID ABUTMENT PLASTIC COPING</b>			
Tissue Level Solid Plastic Coping Crown	ø 5.3	10 mm	S-TSAPICc
Tissue Level Solid Plastic Coping Bridge	ø 5.3	10 mm	S-TSAPICb
<b>TISSUE LEVEL SOLID ABUTMENT IMPRESSION CAP</b>			
Tissue Level Solid Abutment Impression Cap	ø 10.2		S-TSAICap
<b>TISSUE LEVEL SOLID ABUTMENT LAB ANALOG</b>			
Tissue Level Solid Abutment Lab Analog	ø 4.8	14 mm	S-TSALA40
Tissue Level Solid Abutment Lab Analog	ø 4.8	15.5 mm	S-TSALA55
Tissue Level Solid Abutment Lab Analog	ø 4.8	17 mm	S-TSALA70
Tissue Level Solid Abutment Large Lab Analog	ø 6.0	14 mm	S-TSALLA40
Tissue Level Solid Abutment Large Lab Analog	ø 6.0	15.5 mm	S-TSALLA55
Tissue Level Solid Abutment Large Lab Analog	ø 6.0	17 mm	S-TSALLA70
<b>TISSUE LEVEL OCTA SYSTEM</b>			
<b>TISSUE LEVEL OCTA ABUTMENT</b>			
Tissue Level Octa Abutment	ø 3.5		S-TOA
<b>OCTA ABUTMENT HEALING CAP</b>			
Octa Abutment Healing Cap	ø 4.9	4 mm	S-OAHCap
Octa Abutment Healing Cap Large	ø 4.9	6 mm	S-OAHCAP-L
<b>TISSUE LEVEL OCTA SYSTEM PLASTIC CYLINDER</b>			
Octa System Plastic Cylinder Octa	ø 5		S-OAPICy
<b>TISSUE LEVEL OCTA ABUTMENT PLASTIC CYLINDER NON-OCTA</b>			
Tissue Level Octa Abutment Plastic Cylinder Non-Octa	ø 5		S-OAPICyN
<b>TISSUE LEVEL OCTA SYSTEM IMPRESSION COPING PICK UP</b>			
Tissue Level Octa System Impression Coping Pick Up Octa	ø 5.5		S-ICOAP
Tissue Level Octa System Impression Coping Pick Up Non-Octa	ø 5.5		S-ICOAPN
<b>TISSUE LEVEL OCTA SYSTEM IMPRESSION COPING TRANSFER</b>			
Tissue Level Octa System Impression Coping Transfer Octa	ø 5.0		S-ICOAT
Tissue Level Octa System Impression Coping Transfer Non-Octa	ø 5.0		S-ICOATN

<b>PRODUCT</b>	<b>PLATFORM</b>	<b>HEIGHT</b>	<b>CODE</b>
<b>TISSUE LEVEL OCTA ABUTMENT LAB ANALOG</b>			
Tissue Level Octa Abutment Lab Analog	ø 4.8		S-OALA
<b>TISSUE LEVEL O-RING SYSTEM</b>			
<b>TISSUE LEVEL O-RING ABUTMENTS</b>			
Tissue Level O-Ring Abutment	ø 3.5	0.5 mm	S-TORA05
Tissue Level O-Ring Abutment	ø 3.5	2 mm	S-TORA20
Tissue Level O-Ring Abutment	ø 3.5	4 mm	S-TORA40
<b>TISSUE LEVEL O-RING LAB ANALOG</b>			
Tissue Level O-Ring Lab Analog	ø 4.1		S-ORLA
<b>TISSUE LEVEL LOCATOR SYSTEM</b>			
<b>TISSUE LEVEL LOCATOR ABUTMENTS</b>			
Tissue Level Locator Abutment	ø 4.9	1.5 mm	S-TLA15
Tissue Level Locator Abutment	ø 4.9	2.5 mm	S-TLA25
Tissue Level Locator Abutment	ø 4.9	3.5 mm	S-TLA35
Tissue Level Locator Abutment	ø 4.9	4.5 mm	S-TLA45
Tissue Level Locator Abutment	ø 4.9	5.5 mm	S-TLA55
Tissue Level Locator Abutment	ø 4.9	6.5 mm	S-TLA65
<b>LOCATOR MATRIX</b>			
Locator Matrix (includes processing cap)	ø 5.5mm		S-LMHt
Locator Cap (includes 4 pcs caps)	ø 4.86 mm		S-LRI
<b>TISSUE LEVEL LOCATOR LAB ANALOG</b>			
Tissue Level Locator Lab Analog			S-TLLA
<b>TISSUE LEVEL SUPERB OVERDENTURE SYSTEM</b>			
<b>TISSUE LEVEL SUPERB OVERDENTURE ABUTMENTS</b>			
Tissue Level Superb Overdenture Abutment	ø 3.5	0.5 mm	S-TEBA05
Tissue Level Superb Overdenture Abutment	ø 3.5	2 mm	S-TEBA20
Tissue Level Superb Overdenture Abutment	ø 3.5	4 mm	S-TEBA40
<b>SUPERB OVERDENTURE CAP</b>			
Superb Overdenture Cap	ø 5.5 mm	3.5 mm	S-EBOCAP
<b>SUPERB OVERDENTURE LAB ANALOG</b>			
Superb Overdenture Analog	ø 4.1mm	16.6 mm	S-EBLA
<b>ZYGOMA IMPLANT</b>			
<b>ABUTMENTS</b>			
Zygoma Abutment	3 mm		S-ZCA4330
Zygoma Abutment	5 mm		S-ZCA4350

<b>PRODUCT</b>	<b>PLATFORM</b>	<b>HEIGHT</b>	<b>CODE</b>
<b>ZYGOMA ANGLED ABUTMENTS</b>			
Zygoma Angled Abutment		3 mm	S-ZAA43203
Zygoma Angled Abutment		5 mm	S-ZAA43205
<b>CAPS</b>			
Zygoma Healing Cap		3 mm	S-ZHA4330
Zygoma Healing Cap		5 mm	S-ZHA4350
<b>COPING</b>			
Zygoma Impression Coping Pick Up			S-ZICP
<b>PLASTIC CYLINDER</b>			
Zygoma Plastic Cylinder			S-ZPICYN
<b>LAB ANALOG</b>			
Zygoma Lab Analog			S-ZLA
<b>DRIVERS</b>			
Implant Drivers (For ratchet and contra angle)			
Bone Level Mini Fixture Driver Machine Short			S-BMFDMS
Bone Level Mini Fixture Driver Machine Long			S-BMFDML
Bone Level Mini Fixture Driver Ratchet Short			S-BMFDRS
Bone Level Mini Fixture Driver Ratchet Long			S-BMFDRL
Bone Level Fixture Driver Machine Short			S-BFDMS
Bone Level Fixture Driver Machine Long			S-BFDML
Bone Level Fixture Driver Ratchet Short			S-BFDRS
Bone Level Fixture Driver Ratchet Long			S-BFDRL
Tissue Level Fixture Driver Ratchet Short			S-TFDMS
Tissue Level Fixture Driver Ratchet Long			S-TFDRL
Tissue Level Fixture Driver Machine Short			S-TFDMS
Tissue Level Fixture Driver Machine Long			S-TFDML
Hex driver (extra short, short, long, extra long)			S-HDRS / S-HDRL
Hex Driver Contra Angle (extra short, short, long, extra long)			S-HDRS / S-HDRL
Superb Overdenture Driver			S-EBRD
O-ring Driver			S-ORADS/L
Octa Driver			S-OADS/L
Solid Driver			S-SADL
ORRO Insert Adaptor			S-OIA
Bar Matrix Insert Adaptor			S-BMIA

PRODUCT	PLATFORM	HEIGHT	CODE
<b>SCREWS</b>			
<b>CLOSING SCREW</b>			
<b>BONE LEVEL</b>			
Bone Level Mini Closing Screw	ø 2.8	4.5 mm	S-BMCS
Bone Level Closing Screw	ø 3.36	6.5 mm	S-BCS
<b>TISSUE LEVEL</b>			
Tissue Level Closing-Screw Small	ø 3.5	5.35 mm	S-TCSS
Tissue Level Closing-Screw Large	ø 5.2	6.4mm	S-TCSL
Tissue Level Closing-Screw Wide	ø 6.4	6.4 mm	S-TCSW
<b>ABUTMENT SCREWS</b>			
<b>BONE LEVEL</b>			
Bone Level Mini Abutment Screw	ø 2.1	7.64 mm	S-BMAS
Bone Level Abutment Screw	ø 2.3	9.32 mm	S-BAS
Bone Level Multi-Unit Abutment Screw	ø 1.96	9 mm	S-BMUAS
<b>TISSUE LEVEL</b>			
Tissue Level Abutment Screw	ø 2.5	7.1	S-TAS
<b>COMMON SCREWS</b>			
Octa Abutment Screw	ø 2.5	4.5 mm	S-OAS
Multi-Unit Abutment Healing Screw	ø 2.2	4 mm	S-MUAHS



